

new line



a new line of steam sterilizers

Newmed



■ multifunction display
■ controlled by micro-processor

■ double and patented
■ motor-operated closure

■ stainless steel chamber
■ without any welding

■ aluminum trays

■ 2 vacuum pumps

■ built - in printer and
■ usb connection

■ inner heater

■ 16A schuko 230
■ VAC single phase plug


■ osmosis system

complete cabinet




upper sterilization cabinet


134°C class B




121°C class B



134°C prion



134°C class S



an **extra-dry** system can be activated for each sterilization cycle

sterilization temperature:	134°C
sterilization pressure:	2.1 bar
sterilization time:	5'
drying time:	15'
cycle time*:	1h 15'

sterilization temperature:	121°C
sterilization pressure:	1.1 bar
sterilization time:	18'
drying time:	15'
cycle time*:	1h 25'

sterilization temperature:	134°C
sterilization pressure:	2.1 bar
sterilization time:	20'
drying time:	15'
cycle time*:	1h 35'

sterilization temperature:	134°C
sterilization pressure:	2.1 bar
sterilization time:	5'
drying time:	15'
cycle time*:	1h

helix/b&d test

sterilization temperature:	134°C
sterilization pressure:	2.1 bar
sterilization time:	3'30"
cycle time*:	50'

vacuum test

cycle time:	17'
-------------	-----

* Approx. cycle length calculated with warm chamber at full load

comparison table

tiche



■ CLASS B ref. UNI EN 13060	
■ inner heater	
■ overall dimensions/weight: (complete cabinet)	
height:	820 mm (1600mm)
width:	646 mm (680mm)
depth:	783 mm (783mm)
weight:	160kg (205kg)*
■ chamber capacity:	
volume:	59.9l
diameter/depth:	355 mm/605 mm
unwrapped material:	15 kg
wrapped material:	7 kg
■ average water consumption:	3.5l per cycle **
■ power connection:	16A 230VAC 50Hz single phase
	schuko plug

* Approx. value

** Average consumption depending on load and cycle types

standard accessories /options



optional: cabinets for autoclave



■ external labeller



■ different size of sterilization containers:

- 30x30x60 *
- 30x30x30*
- 28.5x46.5x105**
- 28.5x46.5x135**
- 28.5x46.5x150**

optional: label printer and sterilization containers



■ aluminium trays



■ reverse osmosis

standard accessories

Professional Sterilization Systems

Newmed was created to satisfy the needs of customers who demand safety, reliability and simplicity all in one sterilization equipment. And thanks to experience, technology and research, **Newmed** has been able to achieve this goal .



Newmed

Due to a continuous process of development, **Newmed** reserves the right to amend the specification of products and services without notice.



Newmed srl

Via Lenin, 79/A - 42020 Quattro Castella - Reggio Emilia - ITALY
Tel. +39 0522 875166 Fax. +39 0522 243096
info@newmedsrl.it
www.newmedsrl.it