CML Brave Laser

A powerful Laser for multidisciplinary treatments

High performance and high efficiency for dentistry, dermatological surgery and antalgic therapy





RED LASER: 15W - 810 nm

GREEN LASER: 520 nm 2mW



- Periodontology
- Endodontics
- Conservative dentistry
- Prosthesis
- Oral medicine
- Dental aesthetic
- Laser therapy











CML Brave Laser

It is a 15W - 810nm powerful diode Laser suitable for multidisciplinary therapies in dentistry.

The cutting precision, hemostasis, decontamination and biostimulation make diode laser the most versatile instrument in dentistry.

The CML Brave Laser in dental bleaching offers considerable advantages in terms of speed and effectiveness. In the modern dentist practices it becomes a leading tool in the different applications like:

Dental / Aesthetic applications

SURGERY ON SOFT TISSUES

- Bloodless incisions Excision of neoformations
- · Gingivectomy · Frenulectomy, Frenectomy
- Gingivoplasty Abscesses curettage
- Biopsies Haemostasis/coagulation
- Operculectomy Leukoplakia
- · Mucocele/Ranula · Clinical crown lengthening
- Guiding plans and clinical crown exposure to ease eruption
- Removal of Granulation Tissue
- Creation of OVATE PONTICS
- Post-extraction Alveolar Decontamination

CONSERVATIVE DENTISTRY

- Cavity decontamination
- Re-mineralisation of the tooth enamel
- Desensitization of the tooth necks
- Sealing of the dentinal tubules
- Pulp coagulation Vitrification

ENDODONTICS

- Decontamination of the Root Canal
- and periapical lesions resolution
- Pulpotomy

PERIODONTICS

- Decontamination of teeth sulcus
- Periodontal pockets and bifurcations
- Periimplantitis
- Anti-inammatory action against gingivitis

PROSTHESIS

- No retraction cords necessary for the preparation of the gingival sulcus to create the cast
- Sulcus extension

ORAL MEDICINE

- Pre/post-surgery and Pre-implantation bio-stimulation • It treats: aphthous ulceration, Herpes labialis • Lichen planus
- Mycosis Leukoplakia Epulis

DENTAL AESTHETIC

- Tooth discolouration

LASER THERAPY

- Photo-chemical and photo-biological effect in tissue cells • Painless, non-invasive
- Arthro-rheumatic disorders, sports traumatology, rehabilitative therapy, specialist therapies, etc.
- Analgesic rehabilitation
- Acupuncture preparation
- Channelling of active ingredients





PERIODONTICS



ENDODONTICS



CONSERVATIVE DENTISTRY



DENTAL RESTORATIONS



ORAL MEDICINE



DENTAL AESTHETIC

ACCESSORIES AVAILABLE:

- Handpiece for TIP with fiber of 2m
- Arcade TIP for bleaching
- Defocused 45° TIP for dermal treatments and biostimulation
- Fiber Tips of 200 µm straight
- Fiber Tips of 300 µm straight
- Fiber Tips of 400 µm straight
- Optical Fibers of 200µmx2m single/double SMA connector
- Optical Fibers of 400µmx2m single/double SMA connector
- Footswitch
- Wireless Footswitch
- Laser protective glasses
- Plane Wave Tip for Biostimulation
- TIP for Dermal therapy
- "CONTACT" handpiece for surgery
- No "CONTACT" handpiece for bleaching
- Wireless Footswitch
- Power Supply for Lasers
- Suitcase Battery for lasers
- Battery charger
- Remote interlock
- Fiber stripper
- Fiber cutter



TIP HANDPIECE



FIBRE TIPS

INTERCHANGEABLE

BIO-STIMULATION TIP



BLEACHING TIP



WAVE TIP

TECHNICAL FEATURES Max. power on fibre 15W continuos Pulse frequency 0-10 kHz Weight 2.5 kg Wavelength 810 nm & 520 nm Operation mode Display 7" - Colour - Gloves Touch Continuous, Pulsed and Laser Diode AlGaAs, CW (Class IV) High frequency pulsed Laser type **Battery duration** 2 hours at max. power 22 x 24 x 10 cm Certification Pulse Ω 2-1 ms Dimensions CE DOSS

