

# 17 IN 1

## Hydra Facial Machine

 BV LASER

-  Bubble Cleaning
-  Ultrasonic
-  Bipolar HR
-  Micro Current (V Face)
-  Four-pole Roller
-  Spot Clears
-  Color Light Masks
-  Five Poles Colorful Light Microelectronics
-  Cold Hammer
-  Skin Scrubber
-  High Pressure Spray
-  Electrotherapy
-  Plasma
-  Ultrasonic Tip
-  Electric Microneedle
-  Oxygen Revitalizing Bubble
-  Diamond Microdermabrasion



# 8+9 TREATMENT HANDLES

## Standard Handle



Bubble Cleaning



Mico Current (V Face)



Five Poles Colorful Light  
Micro-electronics



Skin Scrubber



High Pressure Spray



Ultrasonic



Bipolar HR



Cold Hammer

## Optional Handle



Oxygen Revitalizing  
Bubble



Spot Clears



Plasma



Diamond  
Microdermabrasion



Electric Microneedle



Ultrasonic Tip



Color Light Masks



Electrotherapy



Four-pole Roller

# ONE MACHINE, MULTIPLE FUNCTIONS



Deep Cleanse Skin



Skin Rejuvenation



Improve Complexion



Reduce Fine Lines



Moisturize Skin



Skin Lifting



Shrink Pores



Deep Cleaning



Blackheads Removing



Skin Firming



Skin Tightening



Skin Care

## FEATURES:


8+9 multifunctional skincare system water oxygen meter can satisfy the customer's different treatment requirements.

Suitable for all skin types, it offers non-invasive, non-surgical treatment options.

4 bottles (A, B, C, D), more convenient for save your time to change bottles constantly.

10.4 large screen, user-friendly interface, easy to operate.

Non-invasive, painless and bloodless therapy.



# TECHNICAL PARAMETERS



Vacuum: >100Kpa

Technology: Hydro dermabrasion, photon light

Max output: 250VA

Treatment handle: 8 (standard) + 9 (optional)

Screen: 10.4" Touch LCD

Voltage: 100-240VAC, 50Hz/60Hz

