



Hubble II Dental Camera Recording System

High-resolution Recording,
Precise Diagnosis & Treatment



Superior Light Source System

Clear Vision, Precise Presentation

- Multiple color temperature adjustments with ultra-high color rendering, accurately restoring tooth colors
- Smart interactive design with touch control + remote control, <0.3 s response, Easy operation in a sterile environment

Details in Full, Worry-free Recording

- HD Image Capture with 1080P Full HD Lens, Every Detail Captured
- Smart Autofocus Technology with 0.1mm Precision, Fast and Accurate Positioning
- Personalized Recording Modes with Multi-Mode Color Display, 40+ Hours of Continuous Recording



HD Large Screen, Real-Time Sharing

- 21.5-inch HD Display
- <0.05s Latency
- Real-time Cross-Department Consultation

Safe, Flexible, and Widely Applicable


- Silent Mobile System
- Stable Load-Bearing Structure
- Multi-Angle Adjustment



Integrated Teaching, Diagnosis, and Communication

- Teaching Demonstration
- Case Management
- Treatment Communication

Specification		
Camera	Resolution	1080P
	Zoom	30 X Optical Zoom
	Auto Focus Range	450mm-1000mm
Dental light	Power input	AC 24V
	Maximum Power	50W
	Illumination	8000lux~40000lux
	color Temperature	White Light:5500k+10% Mixed Light: 4600 k~5200 k, 9 levels adjustable Yellow Light: 3500 k+10%
The Unit	Power input	220V
	Maximum Power	80W



Guangzhou | Los Angeles
Melbourne | Milan

HEADQUARTER

Guangzhou Ajax Medical Equipment Co., Ltd.
Address: Building No.2, Dagang Industrial Zone, Shilou Town, Panyu District, Guangzhou, P.R.China

www.ajaxdent.com
info@ajaxdent.com
+86 020 84847938

+86 18926102494
Ajax Dental
Ajax Dental